



IVTScrip™ mRNA-Human NRP1, (Cap 1, Pseudo-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK5592MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK5592MR is a type of mRNA having 30 nt poly(A) tail and modified with Cap 1 & Pseudo-UTP. It encodes the NRP1 protein. This product can be used in -related researches.

Specifications

Modified bases	Pseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_003873
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

SPECIFICATIONS

Modified bases	Pseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_003873
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

GENE INFORMATION

Alternative Names	BDCA4; NP1
Description	This gene is a protein coding gene. Exhibits GTPase activator activity and vascular endothelial growth factor binding activity. Involved in several processes, including cellular response to hepatocyte growth factor stimulus; positive regulation of cellular component organization; and transmembrane receptor protein tyrosine kinase signaling pathway. Localizes to cytosol; focal adhesion; and plasma membrane.