



## IVTScrip™ mRNA-Human LIN28B, (Cap 1, Pseudo-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK5611MR

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

#### Product overview

This product GTTS-WK5611MR is a type of mRNA having 30 nt poly(A) tail and modified with Cap 1 & Pseudo-UTP. It encodes the LIN28B protein. This product can be used in -related researches.

#### Specifications

<b>Modified bases</b>	Pseudo-UTP
<b>5' Cap</b>	Cap 1
<b>Species</b>	Human
<b>RefSeq</b>	NM_001004317
<b>Applications</b>	Gene therapy research
<b>Format</b>	Powder
<b>Quantity</b>	100 µg
<b>Purification</b>	oligo-dT affinity purification

### SPECIFICATIONS

<b>Modified bases</b>	Pseudo-UTP
<b>5' Cap</b>	Cap 1
<b>Species</b>	Human
<b>RefSeq</b>	NM_001004317
<b>Applications</b>	Gene therapy research
<b>Format</b>	Powder
<b>Quantity</b>	100 µg
<b>Purification</b>	oligo-dT affinity purification

## GENE INFORMATION

<b>Alternative Names</b>	FLJ16517; lin-28B; protein lin-28 homolog B; Lin-28.2; CSDD2
<b>Description</b>	This gene is a protein coding gene. Exhibits RNA binding activity and sequence-specific double-stranded DNA binding activity. Involved in RNA metabolic process and regulation of RNA metabolic process. Localizes to cytosol; nucleolus; and nucleoplasm.