



IVTScrip™ mRNA-Human CD38, (Cap 1, Pseudo-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK5647MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK5647MR is a type of mRNA having 30 nt poly(A) tail and modified with Cap 1 & Pseudo-UTP. It encodes the CD38 protein. This product can be used in -related researches.

Specifications

Modified bases	Pseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001775
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

SPECIFICATIONS

Modified bases	Pseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001775
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

GENE INFORMATION

Alternative Names	ADPRC1; ADPRC 1
Description	This gene is a protein coding gene. Predicted to have NAD ⁺ nucleosidase activity and phosphorus -oxygen lyase activity. Involved in several processes, including B cell receptor signaling pathway; positive regulation of B cell proliferation; and regulation of transcription, DNA-templated. Localizes to extracellular exosome. Implicated in type 1 diabetes mellitus and type 2 diabetes mellitus.