



IVTScrip™ mRNA-Human CD200, (Cap 1, N1-Methylpseudo-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK5709MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK5709MR is a type of mRNA having 30 nt poly(A) tail and modified with Cap 1 & N1-Methylpseudo-UTP. It encodes the CD200 protein. This product can be used in -related researches.

Specifications

Modified bases	N1-Methylpseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_138806.3
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

SPECIFICATIONS

Modified bases	N1-Methylpseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_138806.3
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

GENE INFORMATION

Alternative Names	MRC; MOX1; OX-2; MOX2; CD200
Description	This gene is a protein coding gene. Predicted to have protein binding activity involved in heterotypic cell-cell adhesion. Involved in several processes, including negative regulation of matrix metalloproteinase secretion; positive regulation of arginase activity; and regulation of DNA-binding transcription factor activity. Localizes to axon; cell surface; and neuronal cell body.