



IVTScrip™ mRNA-Human beta-2-microglobulin, (Cap 1, N1-Methylpseudo-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK5793MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK5793MR is a type of mRNA having 30 nt poly(A) tail and modified with Cap 1 & N1-Methylpseudo-UTP. It encodes the beta-2-microglobulin protein. This product can be used in -related researches.

Specifications

Modified bases	N1-Methylpseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_004048
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

SPECIFICATIONS

Modified bases	N1-Methylpseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_004048
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

GENE INFORMATION

Alternative Names	IMD43
Description	<p>This gene is a protein coding gene. Exhibits protein homodimerization activity. Involved in several processes, including amyloid fibril formation; regulation of receptor binding activity; and response to bacterium. Localizes to several cellular components, including Golgi apparatus; HFE-transferrin receptor complex; and MHC class I peptide loading complex. Implicated in arthritis; familial visceral amyloidosis; immunodeficiency 43; and inflammatory bowel disease. Biomarker of several diseases, including artery disease (multiple); encephalomyelitis; measles; obesity; and type 2 diabetes mellitus.</p>