



IVTScrip™ mRNA-Human HLA-ABC, (Cap 1, N1-Methylpseudo-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK5801MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK5801MR is a type of mRNA having 30 nt poly(A) tail and modified with Cap 1 & N1-Methylpseudo-UTP. It encodes the HLA-ABC protein. This product can be used in -related researches.

Specifications

Modified bases	N1-Methylpseudo-UTP
5' Cap	Cap 1
Species	Human
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

SPECIFICATIONS

Modified bases	N1-Methylpseudo-UTP
5' Cap	Cap 1
Species	Human
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

GENE INFORMATION

Alternative Names	HLA-A; HLA-B; HLA-C
Description	This gene is a protein coding gene. The human leukocyte antigen (HLA) system or complex is a group of related proteins that are encoded by the major histocompatibility complex (MHC) gene complex in humans. HLAs corresponding to MHC class I (A, B, and C), all of which are the HLA Class1 group, present peptides from inside the cell.