



IVTScrip™ mRNA-Human PBLD, (Cap 1, N1-Methylpseudo-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK5857MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK5857MR is a type of mRNA having 30 nt poly(A) tail and modified with Cap 1 & N1-Methylpseudo-UTP. It encodes the PBLD protein. This product can be used in -related researches.

Specifications

Modified bases	N1-Methylpseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001033083.1
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

SPECIFICATIONS

Modified bases	N1-Methylpseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001033083.1
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

GENE INFORMATION

Alternative Names	MAWBP; FLJ14767
Description	This gene is a protein coding gene. Exhibits identical protein binding activity. Involved in several processes, including maintenance of gastrointestinal epithelium; negative regulation of epithelial to mesenchymal transition; and negative regulation of transmembrane receptor protein serine/threonine kinase signaling pathway. Localizes to cytoplasm.