



IVTScrip™ mRNA-Human SMAD3, (Cap 1, N1-Methylpseudo-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK5863MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK5863MR is a type of mRNA having 30 nt poly(A) tail and modified with Cap 1 & N1-Methylpseudo-UTP. It encodes the SMAD3 protein. This product can be used in -related researches.

Specifications

Modified bases	N1-Methylpseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001145102
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

SPECIFICATIONS

Modified bases	N1-Methylpseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001145102
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

GENE INFORMATION

Alternative Names	MADH3; HsT17436; LDS1C; hMAD-3; MAD
Description	This gene is a protein coding gene. Exhibits several functions, including SMAD binding activity; enzyme binding activity; and transcription factor binding activity. Involved in several processes, including cell surface receptor signaling pathway; negative regulation of cell cycle; and positive regulation of nitrogen compound metabolic process. Localizes to cytosol; nucleus; and receptor complex. Implicated in Loey's-Dietz syndrome 3; Lynch syndrome; and pancreatic cancer. Biomarker of cervical squamous cell carcinoma; endometrial cancer; and prostate cancer.