



IVTScrip™ mRNA-Human CXCR6, (Cap 1, N1-Methylpseudo-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK5873MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK5873MR is a type of mRNA having 30 nt poly(A) tail and modified with Cap 1 & N1-Methylpseudo-UTP. It encodes the CXCR6 protein. This product can be used in -related researches.

Specifications

Modified bases	N1-Methylpseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_006564
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

SPECIFICATIONS

Modified bases	N1-Methylpseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_006564
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

GENE INFORMATION

Alternative Names	STRL33; TYMSTR; BONZO; CDw186
Description	This gene is a protein coding gene. Predicted to have C-C chemokine binding activity and C-C chemokine receptor activity. Predicted to be involved in several processes, including calcium-mediated signaling; cell chemotaxis; and positive regulation of cytosolic calcium ion concentration. Predicted to localize to external side of plasma membrane. Biomarker of systemic scleroderma.