



IVTScrip™ mRNA-Human CD292, (Cap 1, 5-Methoxy-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK6072MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK6072MR is a type of mRNA having 30 nt poly(A) tail and modified with Cap 1 & 5-Methoxy-UTP. It encodes the CD292 protein. This product can be used in -related researches.

Specifications

Modified bases	5-Methoxy-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_0022112
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

SPECIFICATIONS

Modified bases	5-Methoxy-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_0022112
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

GENE INFORMATION

Alternative Names	ACVRLK3; BMPR-1A
Description	This gene is a protein coding gene. Exhibits several functions, including ATP binding activity; SMA D binding activity; and protein homodimerization activity. Involved in several processes, including p ositive regulation of ossification; positive regulation of transmembrane receptor protein serine/thre onine kinase signaling pathway; and positive regulation of vascular associated smooth muscle cell proliferation. Localizes to external side of plasma membrane. Colocalizes with caveola. Implicated in hereditary mixed polyposis syndrome 2; intestinal disease; and juvenile polyposis syndrome. Bi omarker of pulmonary hypertension.