




IVTScrip™ mRNA-Human CD123, (Cap 1, 5-Methoxy-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK6278MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK6278MR is a type of mRNA having 30 nt poly(A) tail and modified with Cap 1 & 5-Methoxy-UTP. It encodes the CD123 protein. This product can be used in -related researches.

Specifications

Modified bases	5-Methoxy-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001267713.2
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

SPECIFICATIONS

Modified bases	5-Methoxy-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001267713.2
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

GENE INFORMATION

Alternative Names	IL3RY; IL3RX; IL3RAY; hIL-3Ra; CD123; IL-3RA
Description	This gene is a protein coding gene. Predicted to have cytokine receptor activity and interleukin-3 binding activity. Predicted to be involved in cytokine-mediated signaling pathway; monocyte differentiation; and regulation of cell growth. Predicted to localize to external side of plasma membrane and receptor complex.