



IVTScrip™ mRNA-Human OCT3/4, (Cap 1, 2-Thio-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK6389MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK6389MR is a type of mRNA having 30 nt poly(A) tail and modified with Cap 1 & 2-Thio-UTP. It encodes the OCT3/4 protein. This product can be used in -related researches.

Specifications

Modified bases	2-Thio-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001285987.1
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

SPECIFICATIONS

Modified bases	2-Thio-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001285987.1
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

GENE INFORMATION

Alternative Names	Oct-4; Oct-3; OTF3; OTF4
Description	<p>This gene is a protein coding gene. Exhibits several functions, including DNA-binding transcription repressor activity, RNA polymerase II-specific; nucleic acid binding activity; and ubiquitin protein ligase binding activity. Involved in several processes, including BMP signaling pathway involved in heart induction; regulation of canonical Wnt signaling pathway; and regulation of transcription, DNA-templated. Localizes to several cellular components, including cytosol; mitochondrion; and nucleoplasm. Implicated in embryonal carcinoma and seminoma. Biomarker of dysgerminoma; embryonal carcinoma; in situ carcinoma; seminoma; and transitional cell carcinoma.</p>