



IVTScrip™ mRNA-Human CD328, (Cap 1, 2-Thio-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK6407MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK6407MR is a type of mRNA having 30 nt poly(A) tail and modified with Cap 1 & 2-Thio-UTP. It encodes the CD328 protein. This product can be used in -related researches.

Specifications

Modified bases	2-Thio-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_0011362198
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

SPECIFICATIONS

Modified bases	2-Thio-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_0011362198
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

GENE INFORMATION

Alternative Names	AIRM1; CD328; SIGLEC-7
Description	This gene is a protein coding gene. CD328 is a protein that in humans is encoded by the SIGLEC7 gene.