




IVTScrip™ mRNA-Human SOX2, (Cap 1, 2-Thio-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK6503MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK6503MR is a type of mRNA having 30 nt poly(A) tail and modified with Cap 1 & 2-Thio-UTP. It encodes the SOX2 protein. This product can be used in -related researches.

Specifications

Modified bases	2-Thio-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_003106
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

SPECIFICATIONS

Modified bases	2-Thio-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_003106
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

GENE INFORMATION

Alternative Names	ANOP3; MGC2413; SRY-box 2; MCOPS3
Description	This gene is a protein coding gene. Exhibits DNA-binding transcription factor activity; miRNA binding activity; and transcription regulatory region sequence-specific DNA binding activity. Involved in several processes, including endodermal cell fate specification; regulation of protein metabolic process; and sensory organ development. Localizes to cytosol and nucleoplasm. Implicated in breast cancer; nasal cavity cancer; and syndromic microphthalmia 3. Biomarker of breast cancer; melanoma; squamous cell carcinoma; tongue squamous cell carcinoma; and transitional cell carcinoma.