



## IVTScrip™ mRNA-Human CD38, (Cap 1, 2-Thio-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK6547MR

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

#### Product overview

This product GTTS-WK6547MR is a type of mRNA having 30 nt poly(A) tail and modified with Cap 1 & 2-Thio-UTP. It encodes the CD38 protein. This product can be used in -related researches.

#### Specifications

<b>Modified bases</b>	2-Thio-UTP
<b>5' Cap</b>	Cap 1
<b>Species</b>	Human
<b>RefSeq</b>	NM_001775
<b>Applications</b>	Gene therapy research
<b>Format</b>	Powder
<b>Quantity</b>	100 µg
<b>Purification</b>	oligo-dT affinity purification

#### SPECIFICATIONS

<b>Modified bases</b>	2-Thio-UTP
<b>5' Cap</b>	Cap 1
<b>Species</b>	Human
<b>RefSeq</b>	NM_001775
<b>Applications</b>	Gene therapy research
<b>Format</b>	Powder
<b>Quantity</b>	100 µg
<b>Purification</b>	oligo-dT affinity purification

## GENE INFORMATION

<b>Alternative Names</b>	ADPRC1; ADPRC 1
<b>Description</b>	This gene is a protein coding gene. Predicted to have NAD <sup>+</sup> nucleosidase activity and phosphorus-oxygen lyase activity. Involved in several processes, including B cell receptor signaling pathway; positive regulation of B cell proliferation; and regulation of transcription, DNA-templated. Localizes to extracellular exosome. Implicated in type 1 diabetes mellitus and type 2 diabetes mellitus.