




IVTScrip™ mRNA-Human TNFR2, (Cap 1, 2-Thio-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK6571MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK6571MR is a type of mRNA having 30 nt poly(A) tail and modified with Cap 1 & 2-Thio-UTP. It encodes the TNFR2 protein. This product can be used in -related researches.

Specifications

Modified bases	2-Thio-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001066
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

SPECIFICATIONS

Modified bases	2-Thio-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001066
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

GENE INFORMATION

Alternative Names	TNF-RII
Description	This gene is a protein coding gene. Exhibits tumor necrosis factor binding activity and ubiquitin protein ligase binding activity. Involved in cellular response to lipopolysaccharide and positive regulation of membrane protein ectodomain proteolysis. Localizes to membrane raft. Implicated in several diseases, including acne; bone disease (multiple); glomerulonephritis (multiple); lung disease (multiple); and pre-eclampsia. Biomarker of several diseases, including autoimmune disease of musculoskeletal system (multiple); female reproductive organ cancer (multiple); glomerulonephritis (multiple); kidney failure (multiple); and lung disease (multiple).