




IVTScrip™ mRNA-Human HAMP, (Cap 1, 5-Methyl-CTP, 30 nt-poly(A))

Cat. No.: GTTS-WK6747MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK6747MR is a type of mRNA having 30 nt poly(A) tail and modified with Cap 1 & 5-Methyl-CTP. It encodes the HAMP protein. This product can be used in -related researches.

Specifications

Modified bases	5-Methyl-CTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_021175
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

SPECIFICATIONS

Modified bases	5-Methyl-CTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_021175
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

GENE INFORMATION

Alternative Names	LEAP1
Description	This gene is a protein coding gene. Exhibits hormone activity and iron ion transmembrane transporter inhibitor activity. Involved in several processes, including defense response to other organism; iron ion homeostasis; and negative regulation of iron ion transmembrane transport. Localizes to extracellular space. Implicated in hemochromatosis; hemochromatosis type 2B; and hepatocellular carcinoma. Biomarker of several diseases, including anemia (multiple); end stage renal disease; liver disease (multiple); myocarditis; and non-alcoholic fatty liver disease.