



IVTScrip™ mRNA-Human CD203c, (Cap 1, 5-Methyl-CTP, 30 nt-poly(A))

Cat. No.: GTTS-WK6873MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK6873MR is a type of mRNA having 30 nt poly(A) tail and modified with Cap 1 & 5-Methyl-CTP. It encodes the CD203c protein. This product can be used in -related researches.

Specifications

Modified bases	5-Methyl-CTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_005021
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

SPECIFICATIONS

Modified bases	5-Methyl-CTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_005021
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

GENE INFORMATION

Alternative Names	CD203c; PDNP3; NPP3
Description	This gene is a protein coding gene. Exhibits metal ion binding activity; nucleoside-triphosphate diphosphatase activity; and phosphodiesterase I activity. Involved in ATP metabolic process; nucleoside triphosphate catabolic process; and pyrimidine nucleotide metabolic process. Localizes to external side of plasma membrane and perinuclear region of cytoplasm.