



IVTScrip™ mRNA-Human CD111, (Cap 1, 5-Methyl-CTP & Pseudo-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK6978MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK6978MR is a type of mRNA having 30 nt poly(A) tail and modified with Cap 1 & 5-Methyl-CTP & Pseudo-UTP. It encodes the CD111 protein. This product can be used in -related researches.

Specifications

Modified bases	5-Methyl-CTP & Pseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_203286
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

SPECIFICATIONS

Modified bases	5-Methyl-CTP & Pseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_203286
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

GENE INFORMATION

Alternative Names	PRR1; CLPED1; HV1S; HVEC; PRR
Description	This gene is a protein coding gene. Exhibits cell adhesion molecule binding activity and identical protein binding activity. Predicted to be involved in several processes, including heterophilic cell-cell adhesion via plasma membrane cell adhesion molecules; homophilic cell adhesion via plasma membrane adhesion molecules; and positive regulation of viral entry into host cell. Localizes to adherens junction. Implicated in cleft lip; cleft lip-palate-ectodermal dysplasia syndrome; cleft palate; and ectodermal dysplasia.