




IVTScrip™ mRNA-Human ACTR3C, (Cap 0, Pseudo-UTP, 120 nt-poly(A))

Cat. No.: GTTS-WK7374MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK7374MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 0 & Pseudo-UTP. It encodes the ACTR3C protein. This product can be used in Regulatory T (Treg) cell-related researches.

Specifications

Modified bases	Pseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_001164458.2
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	Pseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_001164458.2
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Alternative Names	ARP11
Description	ACTR3C (Actin Related Protein 3C) is a Protein Coding gene. Diseases associated with ACTR3C include Microcephaly, Autosomal Dominant and Lung Adenoma. Among its related pathways are Yersinia infection and Shigellosis. Gene Ontology (GO) annotations related to this gene include actin binding. An important paralog of this gene is ACTR3.