




## IVTScrip™ mRNA-Human ADGRF5, (Cap 0, Pseudo-UTP, 120 nt-poly(A))

Cat. No.: GTTS-WK7440MR

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

#### Product overview

This product GTTS-WK7440MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 0 & Pseudo-UTP. It encodes the ADGRF5 protein. This product can be used in MKI67+ progenitor cell-related researches.

#### Specifications

|                |                       |
|----------------|-----------------------|
| Modified bases | Pseudo-UTP            |
| 5' Cap         | Cap 0                 |
| Species        | Human                 |
| RefSeq         | NM_001098518.1        |
| Applications   | Gene therapy research |
| Format         | Powder                |
| Quantity       | 100 µg                |
| Purification   | Chromatography        |

#### SPECIFICATIONS

|                       |                       |
|-----------------------|-----------------------|
| <b>Modified bases</b> | Pseudo-UTP            |
| <b>5' Cap</b>         | Cap 0                 |
| <b>Species</b>        | Human                 |
| <b>RefSeq</b>         | NM_001098518.1        |
| <b>Applications</b>   | Gene therapy research |
| <b>Format</b>         | Powder                |
| <b>Quantity</b>       | 100 µg                |
| <b>Purification</b>   | Chromatography        |

## GENE INFORMATION

|                    |                                                                                                                                                                                                                                                                                                                                                                                                            |
|--------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b> | ADGRF5 (Adhesion G Protein-Coupled Receptor F5) is a Protein Coding gene. Diseases associated with ADGRF5 include Vibratory Urticaria. Among its related pathways are Metabolism of proteins and GPCRs, Other. Gene Ontology (GO) annotations related to this gene include G protein-coupled receptor activity and transmembrane signaling receptor activity. An important paralog of this gene is ADGRF1. |
|--------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|