




IVTScrip™ mRNA-Human AEN, (Cap 0, Pseudo-UTP, 120 nt-poly(A))

Cat. No.: GTTS-WK7476MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK7476MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 0 & Pseudo-UTP. It encodes the AEN protein. This product can be used in Brush cell (Tuft cell)-related researches.

Specifications

Modified bases	Pseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_022767.4
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	Pseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_022767.4
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Description	AEN (Apoptosis Enhancing Nuclease) is a Protein Coding gene. Gene Ontology (GO) annotations related to this gene include nucleic acid binding and exonuclease activity. An important paralog of this gene is ISG20L2.
--------------------	---