



IVTScrip™ mRNA-Human ABCB8, (Cap 0, N1-Methylpseudo-UTP, 120 nt-poly(A))

Cat. No.: GTTS-WK7542MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK7542MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 0 & N1-Methylpseudo-UTP. It codes the ABCB8 protein. This product can be used in Gonadal mitotic phase fetal germ cell-related researches.

Specifications

Modified bases	N1-Methylpseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_001282291.2
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	N1-Methylpseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_001282291.2
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Alternative Names	MABC1; M-ABC1; MITOSUR; EST328128
Description	This nuclear gene encodes a multi-pass membrane protein that is targeted to the mitochondrial inner membrane. The encoded protein is an ATP-dependent transporter that may mediate the passage of organic and inorganic molecules out of the mitochondria. Loss of function of the related gene in mouse results in a disruption of iron homeostasis between the mitochondria and cytosol. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Aug 2013]