



IVTScrip™ mRNA-Human ABHD14A, (Cap 0, N1-Methylpseudo-UTP, 120 nt-poly(A))

Cat. No.: GTTS-WK7563MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK7563MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 0 & N1-Methylpseudo-UTP. It codes the ABHD14A protein. This product can be used in Retinoid acid signaling-responsive fetal germ cell-related researches.

Specifications

Modified bases	N1-Methylpseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_015407.5
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	N1-Methylpseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_015407.5
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Alternative Names	DORZ1
Description	ABHD14A (Abhydrolase Domain Containing 14A) is a Protein Coding gene. Diseases associated with ABHD14A include Chanarin-Dorfman Syndrome and Autosomal Recessive Non-Syndromic Intellectual Disability. Gene Ontology (GO) annotations related to this gene include hydrolase activity. An important paralog of this gene is ABHD14B.