



IVTScrip™ mRNA-Human ADAMTSL1, (Cap 0, N1-Methylpseudo-UTP, 120 nt-poly(A))

Cat. No.: GTTS-WK7708MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK7708MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 0 & N1-Methylpseudo-UTP. It codes the ADAMTSL1 protein. This product can be used in FGFR1 High NME5-epithelial cell-related researches.

Specifications

Modified bases	N1-Methylpseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_001040272.6
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	N1-Methylpseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_001040272.6
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Alternative Names	C9orf94; PUNCTIN; ADAMTSR1; ADAMTSL-1
Description	This gene encodes a secreted protein and member of the ADAMTS (a disintegrin and metalloproteinase with thrombospondin motif) family. This protein lacks the metalloproteinase and disintegrin-like domains, which are typical of the ADAMTS family, but contains other ADAMTS domains, including the thrombospondin type 1 motif. This protein may have important functions in the extracellular matrix. Alternative splicing results in multiple transcript variants encoding distinct proteins. [provided by RefSeq, Jul 2008]