




IVTScrip™ mRNA-Human ABHD14B, (Cap 0, 2-Thio-UTP, 120 nt-poly(A))

Cat. No.: GTTS-WK8164MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK8164MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 0 & 2-Thio-UTP. It encodes the ABHD14B protein. This product can be used in Liver bud hepatic cell-related researches.

Specifications

Modified bases	2-Thio-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_001146314.2
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	2-Thio-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_001146314.2
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Alternative Names	CIB; HEL-S-299
Description	ABHD14B (Abhydrolase Domain Containing 14B) is a Protein Coding gene. Among its related pathways are Cytochrome P450 - arranged by substrate type and Metabolism. Gene Ontology (GO) annotations related to this gene include hydrolase activity. An important paralog of this gene is ABHD14A.