




IVTScrip™ mRNA-Human AGBL4-IT1, (Cap 0, 2-Thio-UTP, 120 nt-poly(A))

Cat. No.: GTTS-WK8393MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK8393MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 0 & 2-Thio-UTP. It encodes the AGBL4-IT1 protein. This product can be used in Meiotic prophase fetal germ cell-related researches.

Specifications

Modified bases	2-Thio-UTP
5' Cap	Cap 0
Species	Human
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	2-Thio-UTP
5' Cap	Cap 0
Species	Human
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Description	AGBL4-IT1 (AGBL4 Intronic Transcript 1) is an RNA Gene, and is affiliated with the lncRNA class.
--------------------	--