



IVTScrip™ mRNA-Human ABALON, (Cap 0, 5-Methyl-CTP, 120 nt-poly(A))

Cat. No.: GTTS-WK8425MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK8425MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 0 & 5-Methyl-CTP. It encodes the ABALON protein. This product can be used in Mitotic arrest phase fetal germ cell-related researches.

Specifications

Modified bases	5-Methyl-CTP
----------------	--------------

5' Cap	Cap 0
--------	-------

Species	Human
---------	-------

Applications	Gene therapy research
--------------	-----------------------

Format	Powder
--------	--------

Quantity	100 µg
----------	--------

Purification	Chromatography
--------------	----------------

SPECIFICATIONS

Modified bases	5-Methyl-CTP
5' Cap	Cap 0
Species	Human
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Alternative Names	INXS
Description	ABALON (Apoptotic BCL2L1-Antisense Long Non-Coding RNA) is an RNA Gene, and is affiliated with the lncRNA class. Diseases associated with ABALON include Melanoma.