



IVTScrip™ mRNA-Human ABI3BP, (Cap 0, 5-Methyl-CTP, 120 nt-poly(A))

Cat. No.: GTTS-WK8476MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK8476MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 0 & 5-Methyl-CTP. It encodes the ABI3BP protein. This product can be used in FGFR1 High NME5-epithelial cell-related researches.

Specifications

Modified bases	5-Methyl-CTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_001349329.2
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	5-Methyl-CTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_001349329.2
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Alternative Names	TARSH; NESHBP
Description	ABI3BP (ABI Family Member 3 Binding Protein) is a Protein Coding gene. Diseases associated with ABI3BP include Agnathia-Otocephaly Complex and Mood Disorder. Gene Ontology (GO) annotations related to this gene include heparin binding and glycosaminoglycan binding. An important paralog of this gene is FNDC1.