



IVTScrip™ mRNA-Human ACTRT3, (Cap 0, 5-Methyl-CTP, 120 nt-poly(A))

Cat. No.: GTTS-WK8575MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK8575MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 0 & 5-Methyl-CTP. It encodes the ACTRT3 protein. This product can be used in Paneth cell-related researches.

Specifications

Modified bases	5-Methyl-CTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_032487.5
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	5-Methyl-CTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_032487.5
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Description	ACTRT3 (Actin Related Protein T3) is a Protein Coding gene. Diseases associated with ACTRT3 include Dyskeratosis Congenita, Autosomal Dominant 1 and Dyskeratosis Congenita Autosomal Dominant. An important paralog of this gene is ACTB.
--------------------	--