



IVTScrip™ mRNA-Human ADGRF5, (Cap 0, 5-Methyl-CTP & Pseudo-UTP, 120 nt-poly(A))

Cat. No.: GTTS-WK8940MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK8940MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 0 & 5-Methyl-CTP & Pseudo-UTP. It encodes the ADGRF5 protein. This product can be used in MKI67+ progenitor cell-related researches.

Specifications

Modified bases	5-Methyl-CTP & Pseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_001098518.1
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	5-Methyl-CTP & Pseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_001098518.1
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Description	ADGRF5 (Adhesion G Protein-Coupled Receptor F5) is a Protein Coding gene. Diseases associated with ADGRF5 include Vibratory Urticaria. Among its related pathways are Metabolism of proteins and GPCRs, Other. Gene Ontology (GO) annotations related to this gene include G protein-coupled receptor activity and transmembrane signaling receptor activity. An important paralog of this gene is ADGRF1.
--------------------	--