



IVTScrip™ mRNA-Human AC092580.4, (Cap 0, Pseudo-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK9095MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK9095MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 0 & Pseudo-UTP. It encodes the AC092580.4 protein. This product can be used in Natural killer T (NKT) cell-related researches.

Specifications

Modified bases	Pseudo-UTP
5' Cap	Cap 0
Species	Human
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	Pseudo-UTP
5' Cap	Cap 0
Species	Human
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Description	ENSG00000287119 is an RNA Gene, and is affiliated with the lncRNA class. Diseases associated with ENSG00000287119 include Leukemia, Acute Myeloid.
--------------------	--