



IVTScrip™ mRNA-Human AFAP1L2, (Cap 0, Pseudo-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK9280MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK9280MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 0 & Pseudo-UTP. It encodes the AFAP1L2 protein. This product can be used in Oligodendrocyte-related researches.

Specifications

Modified bases	Pseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_001001936.3
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	Pseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_001001936.3
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Description	AFAP1L2 (Actin Filament Associated Protein 1 Like 2) is a Protein Coding gene. Diseases associated with AFAP1L2 include Familial Vesicoureteral Reflux and Vesicoureteral Reflux 1. Among its related pathways are Cytoskeletal Signaling. Gene Ontology (GO) annotations related to this gene include SH3 domain binding and protein tyrosine kinase activator activity. An important paralog of this gene is AFAP1.
--------------------	--