




IVTScrip™ mRNA-Human ABCF2, (Cap 0, N1-Methylpseudo-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK9355MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK9355MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 0 & N1-Methylpseudo-UTP. It codes the ABCF2 protein. This product can be used in Natural killer T (NKT) cell-related researches.

Specifications

Modified bases	N1-Methylpseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_007189.3
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	N1-Methylpseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_007189.3
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Alternative Names	ABC28; HUSSY18; HUSSY-18; EST133090
Description	This gene encodes a member of the ATP-binding cassette (ABC) transporter superfamily. ATP-binding cassette proteins transport various molecules across extra- and intracellular membranes. Alterations in this gene may be involved in cancer progression. Related pseudogenes have been identified on chromosomes 3 and 7. [provided by RefSeq, Mar 2019]