



IVTScrip™ mRNA-Human ABHD17B, (Cap 0, N1-Methylpseudo-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK9367MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK9367MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 0 & N1-Methylpseudo-UTP. It codes the ABHD17B protein. This product can be used in Natural killer T (NKT) cell-related researches.

Specifications

Modified bases	N1-Methylpseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_001025780.3
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	N1-Methylpseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_001025780.3
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Alternative Names	CGI-67; C9orf77; FAM108B1
Description	ABHD17B (Abhydrolase Domain Containing 17B, Depalmitoylase) is a Protein Coding gene. Diseases associated with ABHD17B include Acrofacial Dysostosis 1, Nager Type. Gene Ontology (GO) annotations related to this gene include hydrolase activity and serine-type peptidase activity. An important paralog of this gene is ABHD17A.