



IVTScrip™ mRNA-Human ACSF2, (Cap 0, N1-Methylpseudo-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK9444MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK9444MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 0 & N1-Methylpseudo-UTP. It codes the ACSF2 protein. This product can be used in Leydig cell-related researches.

Specifications

Modified bases	N1-Methylpseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_001288968.2
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	N1-Methylpseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_001288968.2
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Alternative Names	ACSMW; AVYV493
Description	ACSF2 (Acyl-CoA Synthetase Family Member 2) is a Protein Coding gene. Among its related pathways are Fatty Acyl-CoA Biosynthesis and Mitochondrial LC-Fatty Acid Beta-Oxidation. Gene Ontology (GO) annotations related to this gene include ligase activity. An important paralog of this gene is ACSS1.