




IVTScrip™ mRNA-Human ACTRT3, (Cap 0, N1-Methylpseudo-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK9475MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK9475MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 0 & N1-Methylpseudo-UTP. It codes the ACTRT3 protein. This product can be used in Paneth cell-related researches.

Specifications

Modified bases	N1-Methylpseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_032487.5
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	N1-Methylpseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_032487.5
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Description	ACTRT3 (Actin Related Protein T3) is a Protein Coding gene. Diseases associated with ACTRT3 include Dyskeratosis Congenita, Autosomal Dominant 1 and Dyskeratosis Congenita Autosomal Dominant. An important paralog of this gene is ACTB.
--------------------	--