




IVTScrip™ mRNA-Human ADGRG3, (Cap 0, N1-Methylpseudo-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK9542MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK9542MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 0 & N1-Methylpseudo-UTP. It codes the ADGRG3 protein. This product can be used in Natural killer T (NKT) cell-related researches.

Specifications

Modified bases	N1-Methylpseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_001308360.2
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	N1-Methylpseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_001308360.2
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Description	ADGRG3 (Adhesion G Protein-Coupled Receptor G3) is a Protein Coding gene. Diseases associated with ADGRG3 include Polymicrogyria, Bilateral Frontoparietal and Epulis. Among its related pathways are Innate Immune System. Gene Ontology (GO) annotations related to this gene include G protein-coupled receptor activity and transmembrane signaling receptor activity. An important paralog of this gene is ADGRG1.
--------------------	---